

Aneuploidy diagnosis in single cells using array CGH

Dr Michelle Fraser*, Dr Marc Goetz, Tina McCarthy and Dr Nicole Hussey


Reproductive Health Science Pty Ltd, Adelaide, Australia

Dr David Cram, Elissa Osborne and Vicki Kravlevski, Monash IVF and Monash Immunology and Stem Cell Laboratories, Melbourne, Australia

Minimum requirements of array CGH technology for PGS

- Capacity to simultaneously analyse the ploidy status of all chromosomes – 15-30% of aneuploidies are not detected using a standard 9 chromosome FISH probe panel
- Rapid so as to fit in with clinical treatment
- Accurate and reliable diagnostic outcome
- Simple to use and interpret
- Low cost so that large numbers can be processed


Array CGH platforms



- BAC (Bacterial Artificial Chromosome) arrays
- Oligonucleotide arrays
- SNP arrays
- Chromosome library arrays
(RHS array)


The RHS Array

- Diagnosis of every chromosome,
- on single cells,
- within 24 hours



- Contains probes for the 22 human autosomes, X, Y, blank, positive & negative controls each replicated 4 times
- Employs 'DOP' on 'DOP' hybridization

Principles of the RHS Array



Individual microdissected human chromosomes



Unknown test cell




Euploid male reference cell

DOP-PCR amplified



Probe bound to array




Any DOP-PCR bias is shared by all components of the test

Array resolution


- The RHS array looks no further than aneuploidy but:
 - the principle could be applied to smaller regions, and;
 - the DOP products can be used in other tests eg diagnosis of single gene disorders as with other WGA methods

Aneuploid Diagnosis: Trisomy 18



Female Trisomy 18 single fibroblast cell (47,XX+18) versus normal male

Aneuploid Diagnosis: Trisomy 21



Male Trisomy 21 single blood lymphocyte (47,XY+21) versus normal male


Intra and inter-assay variability

Repeat hybridisations of single cell DOP-products:

Karyotype	Number of cells	Number of correct karyotypes
46,XY	2	9/9, 7/7
47,XY,+9	1	4/5 (false negative, 9 = 1.23)
47,XY,+13	1	3/3
47,XY,+15	3	8/8, 6/6, 11/11
47,XX,+18	1	7/7
48,XY,+2,+21	1	9/9
47,XX,+21	2	5/5, 31/33*
45,X	1	9/9
47,XXY	1	6/6
47,XYY	2	2/2, 3/3
47,XXX	1	2/2


Intravariability study: 47,XX,+21

(n=31 hybridisations, 1 cell)




Intravariability study: 47,XX,+21


False Negatives: n=2 hybridisations, 46,XX



Do WGA methods differ?




GenomePlex (Sigma)




Repli-G (Qiagen)

pre-heat fragmentation





post-heat fragmentation





RHS DOP

GenomePlex kit results



Repli-G kit results




PGS study

- 8 aneuploid embryos from 3 patients
- 1 biopsied blastomere diagnosed by FISH
 - X, Y, 13, 18, 21 in first round
 - 15, 16, 22 in second round
- Aneuploid embryos succumbed for 18-24 hours
- Zona dissolved using acid tyrodes and 26 blastomeres transferred to PCR tubes
- 18 of 26 (70%) of blastomeres produced a DOP-PCR product

Embryo 5 (7 blastomeres)

FISH cell 5a = 46,XY,+13


aCGH cell 5b = 46,XY



Embryo 5 (7 blastomeres)

FISH cell 5a = 46,XY,+13

aCGH cell 5c = 48,XY,+17,+19



Overview of Embryo 5

Method	Blastomere	Diagnosis
FISH	a	47,XY,+13
CGH	b	46,XY
CGH	c	48,XY,+17,+19
CGH	d	46,XY
CGH	e	48,XY,+19,+20
CGH	f	46,XY
CGH	g	46,XY
CGH	h & i	failed PCRs

Embryo results

Embryo	FISH result	aCGH results
1	45,XX,-15	45,XX,-4
2	46,XY,-15,+16	47,XY,+1 45,XY,-10,-15,+19 48,XXY,+1,-3,+4,-12,+19
3	46,XX,+16	51,XYY,+5,-6,+12,+14,+19,+22 43,XX,-1,-6,-9,+10,-11,-12,-17 1 x failed PCR
4	46,XX,+13,-18	46,XXY,+13,-14,-19 47,XXY,+13,-14,+16,-19 4 x failed PCR
6	44,XY,-21,-22	42,XY,-3,-9,-19,-22 43,XY,-3,-19,-21 1 x failed PCR
7	45,XX,-15	46,XX,-9,+17
8	44,XX,-15,-21	48,XXY,+17 48,XX,+16,+22 47,XX,+17

Embryo results

Embryo	FISH result	aCGH results
1	45,XX,-15	45,XX,-4
2	46,XY,-15,+16	47,XY,+1 45,XY,-10,-15,+19 48,XXY,+1,-3,+4,-12,+19
3	46,XX,+16	51,XYY,+5,-6,+12,+14,+19,+22 43,XX,-1,-6,-9,+10,-11,-12,-17 1 x failed PCR
4	46,XX,+13,-18	46,XXY,+13,-14,-19 47,XXY,+13,-14,+16,-19 4 x failed PCR
6	44,XY,-21,-22	42,XY,-3,-9,-19,-22 43,XY,-3,-19,-21 1 x failed PCR
7	45,XX,-15	46,XX,-9,+17
8	44,XX,-15,-21	48,XXY,+17 48,XX,+16,+22 47,XX,+17

Summary of FISH/array CGH results on aneuploid embryos

- All 8 embryos were mosaic
- Mean concordance of non-affected FISH chromosomes with array CGH was 75% (range 62.5-82.5%)
- Mean concordance of affected FISH chromosomes with array CGH was 17%
- Additional non-FISH chromosome aneuploidies seen in 15/18 blastomeres (83%)
- Complex abnormalities (> 3 chromosomes) seen in 6/18 blastomeres (33%)
- The data is similar to that reported for FISH re-analysis of aneuploid embryos

Strategies to validate array CGH prior to clinical application

- Hybridise DOP-PCR products using metaphase CGH to confirm results – separate the array performance from the PCR
- Analyse surplus frozen IVF embryos donated to research – young donors, not all aneuploid
- Validate using a range of clinical cell types ie first and second polar bodies and blastocyst