

**ESHRE Cryobiology Mtg – Athens, Greece 9/26/09**

---

# **Vitrification of Oocytes: Biological Lessons Learned From Mice, Applied to Women**



Gary D. Smith Ph.D., HCLD

Associate Professor  
Director of Reproductive Sciences  
Program

Departments of OB/GYN, Physiology, and  
Urology



# Utility of Oocyte Cryopreservation

---

- 1) Preserving fertility in cancer patients and/or women undergoing oophorectomy.
- 2) Ethical / moral / social concerns with embryo cryopreservation.
- 3) Salvaging interrupted IVF cycles.
- 4) Oocyte banking in anticipation of reproduction at an advanced maternal age.

<ul style="list-style-type: none"><li>- Slow-rate freezing</li><li>- Vitrification</li></ul>
--

# Cryo-Damage and Subsequent Oocyte Function




Germinal vesicle-intact (immature)

Metaphase II (mature)


Veeck, 1999

**BIG CELLS**  
Lot of responsibilities




- 1) Spindle Formation / Function
- 2) Microfilament Function
- 3) Zona Pellucida

# Recovery / Survival of Vitrified Mouse MII Oocytes




Morphology


Before Vitrification



After Vitrification




# Influence of Vitrification and Time on Metaphase II Spindle Morphology



# Polyscope® View of Spindle Assembly After Vitriification / Warming


---






# Does Vitricification Affect Spindle Morphology (Length)?


## Strict Criteria:




No measurement






Measurement



# Does Vitrification Cause Aberrant Spindle Chromatin Alignment?



Normal	Abnormal Chromatin Alignment
 <p>Normal MII spindle and polar body. The spindle is clearly visible and aligned.</p>	<p>Control (no vitrification)</p> <p><math>1/37 = 2.7\%^a</math></p>
 <p>Abnormal MII spindle and polar body. An asterisk (*) points to an area of aberrant alignment.</p>	<p>Solution Exposure Only</p> <p><math>1/47 = 2.1\%^a</math></p>
 <p>Abnormal MII spindle and polar body. An asterisk (*) points to an area of aberrant alignment.</p>	<p>Vitrification</p> <p><math>3/49 = 6.1\%^a</math></p> <p><math>^aP=0.63</math></p>



# Spindles and Cryopreservation



---

## The Bottom Line:


- 1) All three forces are detrimental
  - Mechanical (Ice crystals)
  - Thermal (Microtubule depolymerization)
  - Chemical (DMSO & PROP)
- 2) Debate on which meiotic stage most susceptible to damage
- 3) Remember microtubules are in dynamic flux
  - polymerization / depolymerization

Today: greatest success is cryopreservation of MII oocytes

# Oocyte Vitrification and Fertilization By Standard Insemination




# Why Do Vitrified Oocytes Not Fertilize?




# Does Vitrification Induce Premature Cortical Granule Release?


---



# Can Fetuin in Cryo-Media Influence Zona Pellucida Modifications?



# Mouse Oocyte Vitrification, ICSI, Fertilization, and Live-Births



# Oocyte Cryopreservation

---

## 1) Slow-rate freezing can work

- 50 to 76% survival

- ~13 to 35% pregnancy rate/transfer

  - (Porcu *et al.*, 2000, *Mol & Cell Endo*)

  - (Boldt *et al.*, 2006, *Reprod Biomed Online*)

  - (Bianchi *et al.*, 2006, *Reprod Biomed Online*)

## 2) Vitrification current area of focus

- two primary factors to consider for success

  - a - vitrification means (container and coolant)

  - b - solutions



# Vitrification Containers

---

## What makes a good vitrification container?

- 1) Low volume of cryo-solution
  - rapid heat transfer
  - fastest cryopreservation / warming
  - no ice crystal formation
- 2) User friendly
- 3) Security

EM Grids - Martino et al., 1996, *Biol Reprod*, Vol 54



Open Pulled Straws (OPS) - Vajta et al., 1998, *Mol Reprod & Dev*, Vol 51

Cryoloops® - Lane et al., 1999, *Fertil & Steril*, Vol 72

Cryotop® - Katayama et al., 2003, *Fertil & Steril*, Vol 80

CryoTip® - Kuwayama and Smith

# Closed Pulled Straw System



- low volume  $\sim 0.4 \mu\text{l}$   
(inner diam. =  $200 \mu\text{m}$ )
- secure and protected  
(heat-sealed at both ends)

# Vitrification / Warming Solutions

---

## Vitrification Solutions

7.5% ethylene glycol  
7.5% DMSO  
20% Protein  
in M199-H

Equilibration  
Solution



15% ethylene glycol  
15% DMSO  
0.5 M sucrose  
20% Protein  
in M199-H

Vitrification  
Solution

## Warming Solutions

### Initial Warming Solution (IWS)

- M199-H
- 1.0 M Sucrose
- 20% Protein

### Dilution Solution (DS)

- M199-H
- 0.5 M Sucrose
- 20% Protein

### Wash Solution (WS)

- M199-H
- 20% Protein

# Objective

Compare human metaphase II oocyte cryopreservation by slow-rate freezing and vitrification in a prospective, randomized, controlled clinical study.





**Joyce Fioravanti, BSc**

2004 Technical Training:

- 1) Oocyte slow-rate freezing
  - Dr. Porcu, M.D.
  - Bologna, Italy
- 2) Oocyte vitrification
  - University of Michigan


# Prospective Randomized Controlled Study

- IRB approval and written informed consent
- Offered to women > 9 oocytes collected in fresh cycles  
Randomized to slow-rate freezing (Porcu) vs. vitrification
- Began January, 2005
- As of Dec 2008: 230 cases cryopreserved, 78 thaws/warmings
- Four hours incubation before intracytoplasmic sperm injection






# Survival and Fertilization




	Age	Thaws/ Warming	Survival	Normal Fertilization	Cleavage
Slow-Rate Freezing	31±1	30	67% <sup>a</sup>	67% <sup>c</sup>	70% <sup>a</sup>
Vitrification	32±1	48	81% <sup>b</sup>	77% <sup>d</sup>	91% <sup>b</sup>

# Day 3 Embryo Development Following Oocyte Cryopreservation





# Prospective Randomized Controlled Study: Pregnancy Rates



# Oocyte Cryopreservation Meta-Analysis

---

	Vitrification		Slow-Rate Freezing	
Age	32	32	34	31
Fertilization Rate	74%	77%	65%	67%
Clinical PR Per Oocytes Thawed/Warmed	4.5%	5.2%	2.3%	1.7%
Clinical PR Per Thawing/Warming	?	38%	?	13%
Clinical PR/ Transfer	46%	38%	21%	21%

Oktay et al., *Fertil & Steril*, 2006

Meta-Analysis

Prospective-Randomized

Meta-Analysis

Prospective-Randomized

# Conclusions

---

- 1) Oocyte cryopreservation by slow-rate freezing or vitrification are experimental procedures to be performed under IRB and informed consent.
- 2) Long-term follow-up of offspring obtained through oocyte cryopreservation is wanting (human and model systems).
- 3) Vitrification is not difficult.
- 4) Vitrification requires practice, can carry a technical signature, and can be “unforgiving” to technical mistakes.

# Acknowledgements

---

## **Program Support**

Timothy Johnson, MD

Max Wicha, MD

## **Physician Support**

Gregory Christman, MD

Senait Fisseha, MD

Eduardo LA Motta, MD, PhD

Dana Ohl, MD

John Randolph, MD

Paulo Serafini, MD, PhD

Yolanda Smith, MD

## **Research Support**

Cristine Silva E Silva

Joyce Fioravanti, BS

Irvine Scientific

Beto Alegretti, BS

Pericles Hassun, PhD

Keiko Hata, MS

Satoshi Kishigami, PhD

Masa Kuwayama, PhD

Marcia Leonard, RN, PNP

Thom Saunders, PhD

Jason Swain, PhD

**National Institutes of Health - NICHD and NCI**



**Thank You**  
[smithgd@umich.edu](mailto:smithgd@umich.edu)